

CreaTec's Bake Unit is used to control the bake-out process of vacuum systems. It contains several independent power lines for the different heating stages, a Eurotherm controller for temperature control via type K thermocouple and a timer for time control.



TYPE	BU-3-230-23A
INPUT VOLTAGE	3 phases 400 V AC
EUROTHERM CONT.	2116
OUTPUT VOLTAGE	6 (or 12) x 230 V AC
MAX. OUTPUT POWER	22 kW
THERMOCOUPLE	type K or as specified
INTERFACE	5V TTL logic signal
ALARM	standard
TIME INTERVAL	0 - 60 h
SIZE	19" x 3 HE x 545 mm