

C6657 DSP Upgrade

This upgrade set allows users of the SBC6711 electronics (USB connection to the PC) to upgrade the DSP to the newer much more powerful C6657 model (Ethernet connection to the PC). The new C6657 DSP allows running the most recent version of the software including more powerful scripting features. The hardware limits of the analogue I/O persist as the complete analogue part of the electronics remains untouched.

The new software version 4.x has several new features like: Automated boot up via integrated DHCP- and TFTP-server, Form and Log info, Scancontrol, second independent Data Recorder, integrated Scripting window and many more. Using the (new) integrated STM simulation (running on the PC) one may use an existing *.dat file to try all software functions before purchasing the DSP upgrade. Just log into your existing account and choose DSP Simulator during the installation of version 4.x.

NOTE



The oldest version of the SBC6711 hardware may not be upgraded, because it does not have the special HD50 socket (LDVS Link) in the backplane.

TYPE	C6657 DSP Upgrade
DIMENSIONS (LxWxH)	220x165x52
CABLE SET	- Power Supply - Ethernet connection to PC - HD50/LDSV Link

