

CreaTec's Shutter Controller is used to control up to 16 channels with 24 UDC signals. Either electrical device or pneumatic actuator (with valve clusters or single valves) are possible. For each channel, the sensor of the correspondent device will be read in real time. Furthermore data is logged in real time and can be exported to an USB memory stick.

Features

- Up to 16 channels
- 24 UDC signals
- Front panel operation (touch screen) and serial interface (RS-232)
- Real time sensor reading and actuator activation
- Data logging and data export to USB storage

TYPE	EPS
INPUT VOLTAGE	100 - 240 VAC
INTERFACE	RS 232
OPERATION	front touch panel
SIZE	19" x 3 HE
POWER CONSUMPTION	300 W