

CreaTec offers custom designed quartz crystal microbalance sensors to measure the flux density. The sensor head is water cooled and the gold plated quartz plates can be easily exchanged. Power supply with two inputs, control unit and oscillator with all necessary cables are included.

## Options

- Linear stroke
- Linear or rotary shutter
- Pneumatic/electric actuation of stroke or shutter
- cVac software control
- Spare quartz crystal

TYPE	QCM
FLANGE SIZE	DN 40 CF   63 or as specified
MOUNTING DIRECTION	as specified
MAX. FREQUENCY	6 MHz
FREQUENCY RESOLUTION	0,1 Hz   0,01Hz
INTERFACE	RS 32   485; ethernet

